

# Determination of molecular formula of a compound using spectroscopy

Abhishek Dasgupta (06MS07)

Indian Institute of Science Education and Research, Kolkata

April 24, 2007

# Introduction

- ▶ My project: Find out molecular formula of **ferric salicylate**
- ▶ I used Job's method for doing this experiment.

# Job's Method

- ▶ Method of continuous variation.
- ▶ Dependent on validity of Beer's Law

$$A = \log I_0/I$$

and

$$A = \epsilon cl$$

# Job's Method

- ▶ Method of continuous variation.
- ▶ Dependent on validity of Beer's Law

$$A = \log I_0/I$$

and

$$A = \epsilon cl$$

- ▶ We vary the concentration of  $\text{Fe}^{3+}$  and salicylic acid so that  $[\text{Fe}^{3+}] + [\text{Salicylic acid}]$  is constant.

# Job's Method

- ▶ Method of continuous variation.
- ▶ Dependent on validity of Beer's Law

$$A = \log I_0/I$$

and

$$A = \epsilon cl$$

- ▶ We vary the concentration of  $\text{Fe}^{3+}$  and salicylic acid so that  $[\text{Fe}^{3+}] + [\text{Salicylic acid}]$  is constant.
- ▶ ... in the ratios 1:9, 2:8, 3:7, ... 9:1.
- ▶ Get spectral scans 400 nm - 650 nm; visible region
- ▶ See peaks; plot absorbance vs composition of solution.
- ▶ Maxima gives us molecular composition!

# Procedure

## What I used

A Spectrophotometer, Ferric Nitrate, Salicylic Acid, Glass beakers (100 mL), Round bottomed flasks (3)

# Procedure

## What I used

A Spectrophotometer, Ferric Nitrate, Salicylic Acid, Glass beakers (100 mL), Round bottomed flasks (3)

## Stocks

- ▶ 0.0025M Sulphuric acid (Solvent). Prepared by standardisation with NaOH, which in turn was standardised against oxalic acid. **Reference spectrum**

# Procedure

## What I used

A Spectrophotometer, Ferric Nitrate, Salicylic Acid, Glass beakers (100 mL), Round bottomed flasks (3)

## Stocks

- ▶ 0.0025M Sulphuric acid (Solvent). Prepared by standardisation with NaOH, which in turn was standardised against oxalic acid. **Reference spectrum**
- ▶ 0.0025M Ferric nitrate, 0.0025M Salicylic acid.



# Procedure

## What I used

A Spectrophotometer, Ferric Nitrate, Salicylic Acid, Glass beakers (100 mL), Round bottomed flasks (3)

## Stocks

- ▶ 0.0025M Sulphuric acid (Solvent). Prepared by standardisation with NaOH, which in turn was standardised against oxalic acid. **Reference spectrum**
- ▶ 0.0025M Ferric nitrate, 0.0025M Salicylic acid.

Stocks showed negligible absorbance.

## Procedure *contd.*

- ▶ Mix stocks  $\rightarrow$  purple colour. Complex formation.


## Procedure *contd.*

- ▶ Mix stocks  $\rightarrow$  purple colour. Complex formation.
- ▶ 9 compositions taken: varying from 10 percent to 90 percent of Ferric nitrate.




- ▶ Spectral scan of visible region.

# Spectra



# Spectra



# Observation


As we see we get maximum absorbance at 5:5 at about 525nm.

# Observation

As we see we get maximum absorbance at 5:5 at about 525nm. To illustrate this we pick 525 nm, and plot absorbance vs composition.

# Observation

As we see we get maximum absorbance at 5:5 at about 525nm. To illustrate this we pick 525 nm, and plot absorbance vs composition.





# Conclusion

From the graph we clearly see that maximal absorption at **1:1** ratio.  
So we have got the molecular formula!

# Conclusion

From the graph we clearly see that maximal absorption at **1:1** ratio.  
So we have got the molecular formula!

## Precautions

- ▶ Proper standardisation

# Conclusion

From the graph we clearly see that maximal absorption at **1:1** ratio.  
So we have got the molecular formula!

## Precautions

- ▶ Proper standardisation
- ▶ Spectral data should be in range where Beer's law is valid; otherwise we absorbance will not have a linear relationship with concentration.

# Conclusion

From the graph we clearly see that maximal absorption at **1:1** ratio.  
So we have got the molecular formula!

## Precautions

- ▶ Proper standardisation
- ▶ Spectral data should be in range where Beer's law is valid; otherwise we absorbance will not have a linear relationship with concentration.

..and so my project ended.

# Other things that could be done

## The stability constant

... is given by

$$K = \frac{[\text{complex}]}{[\text{Fe}^{3+}][\text{Sal}]}$$

We already know  $[\text{Fe}^{3+}]$  and  $[\text{Sal}]$ ; if we could find out the extinction coefficient  $\varepsilon$  for  $[\text{complex}]$ , then we would be done.

# Acknowledgements

I thank Dr. Pradip Bag and Dr. Parna Gupta Bhattacharya, whose guidance greatly helped me in the completion of my project.

I thank our Director, Dr. Sushanta Dattagupta, for giving us encouragement and organising the necessary infrastructure for my project to be possible.

I also thank Mr. Sudhangshu Maity and everyone in the Department of Chemistry for helping me in executing this project.